



CLEAN AIR



CLEAN AIR effectively eliminates excess halogenated anesthetic gases after patient inhalation in the operating room environment.

CLEAN AIR offers an average useful life of 12 to 15 hours.

CLEAN AIR can remove no fewer than 50 grams of pure halogenated anesthetic gases and can be conveniently discarded.

WARNING: NOT FOR USE WITH FLAMMABLE GASES.

Keep from fire and excessive heat.

Submerge in water to reduce flammability hazard if inadvertently used with flammable gases. Container should not be punctured or incinerated.

| ITEM # | DESCRIPTION | UNIT |
|--------|-------------|-------|
| 6000 | CLEAN AIR | 1 Can |