




# ANESORB



**ANESORB effectively eliminates excess halogenated anesthetic gases after patient inhalation in the operating room environment.**

ANESORB offers an average useful life of 12 to 15 hours.

ANESORB can remove no fewer than 50 grams of pure halogenated anesthetic gases and can be conveniently discarded.

**WARNING: NOT FOR USE WITH FLAMMABLE GASES.**

Keep from fire and excessive heat.

Submerge in water to reduce flammability hazard if inadvertently used with flammable gases. Container should not be punctured or incinerated.

ITEM #	DESCRIPTION	UNIT
7000	ANESORB	1 Can